

# Luvis E100 Clamp Mounted LED Exam Light



*Stock Photo - Options May Vary*

## FACTORY NEW

The Luvis E100 LED Exam Light is ideal for medical examinations in smaller spaces. The E100 delivers value through more extended life, lower energy consumption, and will save the cost of hassled maintenance.

The E100 LED Exam Light features three selectable color temperature options (3,500K, 4,000K, and 4,500K). The Luvis E100 Exam Light is available in ceiling, mobile, and wall configurations and comes with a 2-year mfg. Warranty.

Are you looking for a specific LED Surgical Light or additional types of surgical lights? Click any of the following links to view our complete catalog of LED Surgical Lights, Surgical/OR Lights, Minor Surgical/Procedural Lights, or Exam/Diagnostic Lights.

## FEATURES

- Excellent Illumination
- Ergonomic & Friendly Design
- No UV Emission & Cool Light
- Excellent Depth of Illumination
- 3 Selectable Color Temperature Options

## SPECIFICATIONS

- LEDs: 8
- CRI: 95
- Color Temperature
- 3,500K /4,000K /4,500K
- Head Size: 9.84" x 9.68"
- Focus Size: 7" (Circle)/11.8" (Circle)

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free      Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)