

# Luvis E100 LED Exam Light | Wall Mounted



Stock Photo - Options May Vary

## FACTORY NEW

The Luvis E100 LED Exam Light is ideal for medical examinations in smaller spaces. The E100 delivers value through longer life, lower energy consumption, and will save the cost of hassled maintenance. The E100 LED Exam Light features three selectable color temperature options (3,500K, 4,000K, and 4,500K). The Luvis E100 Exam Light can come in ceiling, mobile, and wall configurations and comes with a 2-year mfg.

warranty. Have Questions? Contact us today for questions or available options to fit your specific needs!

## FEATURES

- Excellent Illumination
- Ergonomic & Friendly Design
- No UV Emission & Cool Light
- Excellent Depth of Illumination
- 3 Selectable Color Temperature Options

## SPECIFICATIONS

- LED's: 8
- CRI: 95
- Color Temperature: 3,500K /4,000K /4,500K
- Head Size: 9.84" x 9.68"
- Focus Size: 7" (Circle)/11.8" (Circle)

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)