

MCM Lift Assist 6-Leg/28 Hook I.V. Pole | MCM 253



Stock Photo - Options May Vary

FACTORY NEW

MCM 253 Lift Assist I.V. pole. Mid Central Medical is an American company committed to manufacturing the highest quality medical equipment to serve the demands of America's healthcare industry. In addition, the MCM253 Lift Assist I.V.

Pole (shown) with 28-Hook Top, easily lifts heavy loads by simply turning the crank, eliminating back injuries! Additional MCM 253 features include height adjustability from 56" to 96" and a stable 6-leg base. It's 3" scuff-resistant casters provide smooth maneuverability.

Additionally, it requires no electric motors or hydraulic pumps, reducing maintenance costs and outlasting competitors. It also comes with optional 4-hook and 8-hook tops.

FEATURES

- Easily Lifts Heavy Loads
- 28-Hook Top
- No Electric Motors or Hydraulic - Simply Turn Crank
- Eliminates Back Injuries
- 6-Leg Base

SPECIFICATIONS

- 28-Hook Top
- Manual Turn Crank
- 6-Leg Base
- Height: 56"-96"
- Casters: 3" Easy Maneuverability

**Brochure and images are a general overview of the product. Actual product options may vary.*