

# MCM Infusion Pump Frame I.V. Pole Accessory | MCM 259



Stock Photo - Options May Vary

## FACTORY NEW

The MCM 259 IV Infusion Pump Frame is designed to securely hold up to four infusion pumps. Built with 1" diameter tubular chrome construction, it provides durability and stability in demanding healthcare environments while maintaining a professional appearance.

This IV pole accessory is compatible with most standard setups and requires two MCM 231 Universal Clamps for secure mounting. The chrome finish resists buildup and cleans easily, making maintenance straightforward in busy clinical settings.

Ideal for hospitals, emergency rooms, and outpatient clinics, the MCM 259 helps healthcare professionals organize and support multiple infusion pumps efficiently while saving valuable floor space.

## FEATURES

- Construction: 1" Diameter Tubular Chrome
- Designed to Easily Adapt Multiple Infusion Pumps
- Easily Mounts with two MCM 231 Universal Clamps (required)
- Frame Holds up to 4 Pumps

## SPECIFICATIONS

- Construction: 1" Diameter Tubular Chrome
- Frame Holds up to 4 Pumps
- Frame Measures 3/4" x 27 3/4"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)