

MCM Lift Assist 6-Leg/8 Hook I.V. Pole | MCM 252



Stock Photo - Options May Vary

FACTORY NEW

MCM 252 Lift Assist IV Pole. Mid Central Medical is an American company committed to manufacturing the highest quality medical equipment to serve the demands of America's health care industry. Furthermore, the 252 Lift Assist I.V.

Pole (shown) with 4-Hook Top, easily lifts heavy loads by simply turning the crank, eliminating back injuries! Additionally, the MCM 252 Features: Height is adjustable from 56" - 96", and the 6-leg base offers the stability you need. 3" scuff-resistant casters provide unmatched maneuverability. Moreover, no electric motors or hydraulic pumps to be concerned about or replaced. Saves money initially and will outlast its peers.

FEATURES

- Easily Lifts Heavy Loads
- No Electric Motors or Hydraulic - Simply Turn Crank
- Eliminates Back Injuries
- 8-Hook Top
- 6-Leg Base
- Height Adjustable from 56"-96" High

SPECIFICATIONS

- 8-Hook Top
- Manual Turn Crank
- 6-Leg Base
- Height: 56"-96"
- Casters: 3" Easy Maneuverability

**Brochure and images are a general overview of the product. Actual product options may vary.*