

# MCM Lift Assist 6-Leg/4 Hook I.V. Pole | MCM 251



Stock Photo - Options May Vary

## FACTORY NEW

MCM 251 Lift Assist IV Pole. MCM is a company committed to manufacturing the highest quality medical equipment to serve the demands of America's healthcare industry. Additionally, the MCM 251 Lift Assist I.V.

Pole w/4-Hook Top easily lifts heavy loads by simply turning the crank, eliminating back injuries! Additional MCM 251 Features: Height is adjustable and the 6-leg base offers the stability you need. Furthermore, the scuff-resistant casters on the MCM 251 provide unmatched maneuverability.

No electric motors or hydraulic pumps to be concerned about or replaced. Moreover, the MCM 251 saves money initially and will outlast its peers. Also available with an 8-hook top and a 28-hook top (not shown).

## FEATURES

- Easily Lifts Heavy Loads
- 4-Hook Top
- No Electric Motors or Hydraulic - Simply Turn Crank
- Eliminates Back Injuries
- 6-Leg Base

## SPECIFICATIONS

- 4-Hook Top
- Manual Turn Crank
- 6-Leg Base
- Height: 56"-96"
- Casters: 3"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)