

MCM 12 Qt. Kick Buckets/Sponge Receptacle | MCM 540



Stock Photo - Options May Vary

FACTORY NEW

MCM 540 Kick Buckets/Sponge Receptacle. Looking for dependable and well-designed surgical support equipment? The MCM 540 Kick Bucket/Sponge Receptacle is a high-performing, space-efficient solution trusted by medical professionals across a range of clinical environments.

The 12 qt. capacity handles the rigors of daily use while offering easy handling and simple sanitation. Whether in the OR or procedural room, it delivers reliable functionality where it matters most.

Crafted with a durable stainless steel frame, the MCM 540 is equipped with a full perimeter bumper to help prevent wall and furniture damage during movement. It also features 2" swivel scuff-resistant casters, providing smooth, quiet, and precise maneuverability across different flooring types.

FEATURES

- Stainless steel frame
- Perimeter bumper to protect walls
- Swivel scuff-resistant casters
- Removable Basin

SPECIFICATIONS

- Construction: Stainless Steel
- Dimensions: 14"H x 13 1/2" diameter
- Capacity: 12 Quarts

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)