

# MCM 8.5 Qt. Kick Bucket/Sponge Receptacle | MCM 541



Stock Photo - Options May Vary

## FACTORY NEW

MCM 541 Kick Bucket/Sponge Receptacle - 8.5Qt. Mid Central Medical continues to design and manufacture innovative, high-quality medical equipment to meet the evolving demands of today's healthcare professionals. The MCM 541 Kick Bucket reflects our commitment to durability, functionality, and reliable performance.

We craft every product with precision, using premium materials that support daily clinical use and exceed customer expectations. The MCM 541 Kick Bucket is constructed with a strong stainless steel frame, designed to withstand repeated use in surgical and clinical environments.

Its sleek, polished design resists corrosion and provides a professional appearance suitable for any medical setting.

## FEATURES

- Stainless steel frame
- Perimeter bumper to protect walls
- Swivel scuff-resistant casters
- Removable Basin

## SPECIFICATIONS

- Construction: Stainless Steel
- Dimensions: 8"H x 16-1/4" diameter
- Capacity: 8.5 Quarts

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)