

# MCM Stainless Steel Solution Stand | MCM 1001



Stock Photo - Options May Vary

## FACTORY NEW

MCM 1001 - Stainless Steel Solution Stand. The MCM 1001 from Mid Central Medical is strong and reliable for busy medical settings. As an American company, Mid Central Medical uses quality materials and skilled work to make sure it works well and lasts a long time.

This single basin stand with a shelf is made from strong 16-gauge #304 stainless steel. Additionally, it has a large 8½-quart basin that holds plenty of solutions and supplies. Furthermore, the extra shelf gives space to store tools and equipment.

Moreover, with 3-inch scuff-resistant casters, the stand moves smoothly and quietly on clinic floors, making it easy to use and move around.

## FEATURES

- Durable Stainless Steel Construction
- Stainless Steel Shelf
- Large Basin
- Rolls quietly along on 3" Scuff-Resistant Casters
- Rubber bumpers around the basin stop noise when rolling

## SPECIFICATIONS

- Dimensions: 15" x 18" x 34"
- Construction: 16 Gauge #304 Stainless Steel
- Basin Capacity: 8-1/2 qt

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)