

MCM 20 x 16 Instrument Table w/H Brace | MCM 510



Stock Photo - Options May Vary

FACTORY NEW

The MCM 510 is a heavy-duty stainless steel instrument table designed for medical facilities requiring durable, reliable equipment. Constructed from 16-gauge #304 stainless steel with solid welded construction, the table features square tube legs with professionally welded and polished corners for maximum strength and longevity. The 1.5-inch revealed edge provides superior structural integrity.

The table measures 20 inches long by 16 inches wide by 34 inches high, making it suitable for standard operating room and clinical environments. Standard 3-inch ball bearing casters enable smooth mobility, while the seamless design eliminates open cavities that could harbor contaminants.

Built to withstand demanding healthcare settings, this instrument table combines functional design with corrosion-resistant materials. The unit comes with a limited 2-year warranty and is available in various sizes and configurations to meet different facility requirements.

FEATURES

- Heavy Duty 16 gauge #304 stainless steel with a solid welded construction
- The revealed edge is 1 1/2" to give you the strongest table known
- Square tube legs allow more welding surface which has proven to outlast any round tube leg
- Corners are professionally welded and polished smooth
- 3" ball bearing casters are standard on all tables
- Reinforcing channels welded into larger models for additional strength to meet the highest expectations
- No open cavities to collect unwanted substances
- Cost-effective to meet your budget!

SPECIFICATIONS

- Dimensions: 20"L x 16"W x 34"H
- Construction: Stainless Steel, Solid Weld
- Casters: Ball Bearing
- Caster Diameter: 3"
- Warranty: Limited

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com