

MCM 20 x 16 Instrument Table w/Shelf | MCM 501



Stock Photo - Options May Vary

FACTORY NEW

The MCM 501 is a compact instrument table constructed from 16-gauge #304 stainless steel with solid welded construction. Designed for surgical centers, clinics, and hospitals, it provides a durable solution for facilities with limited space. The sturdy design ensures long-term reliability and easy cleaning.

The table features square tube legs with professionally welded and polished corners, providing superior strength and durability. Standard 3" ball bearing casters enable easy mobility throughout the facility. The design includes no open cavities to prevent collection of unwanted substances, meeting strict hygiene standards.

An integrated shelf provides convenient storage for essential instruments and supplies. The reinforced construction and solid welded design make this table suitable for demanding healthcare environments requiring dependable performance.

FEATURES

- Heavy Duty 16-gauge #304 stainless steel with a solid welded construction
- The revealed edge is 1-1/2" to give you the strongest table known
- Square tube legs allow more welding surface which has proven to outlast any round tube leg
- Corners are professionally welded and polished smooth
- 3" ball bearing casters are standard on all tables
- Reinforcing channels welded into larger models for additional strength to meet the highest expectations
- No open cavities to collect unwanted substances
- Cost-effective to meet your budget!

SPECIFICATIONS

- Dimensions: 20"L x 16"W x 34"H
- Construction: Stainless Steel, Solid Weld
- Casters: Ball Bearing
- Caster Diameter: 3"
- Warranty: Limited

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)