

MCM 33 x 18 Instrument Table w/H Brace | MCM 512



Stock Photo - Options May Vary

FACTORY NEW

The MCM 512 is a heavy-duty stainless steel instrument table designed for medical facilities requiring durable, reliable equipment. Constructed from 16-gauge #304 stainless steel with solid welded construction, the table features a 1 1/2" revealed edge and square tube legs with professionally welded and polished corners for maximum strength and durability.

The table is equipped with 3" ball bearing casters for easy mobility and includes an H-brace for additional structural support. The design eliminates open cavities that could collect unwanted substances, making it easy to maintain and sanitize for clinical environments.

This instrument table is built to withstand the demands of healthcare settings and comes with a limited 2-year warranty.

FEATURES

- Heavy Duty 16 gauge #304 stainless steel with a solid welded construction
- The revealed edge is 1 1/2" to give you the strongest table known
- Square tube legs allow more welding surface which has proven to outlast any round tube leg
- Corners are professionally welded and polished smooth
- 3" ball bearing casters are standard on all tables
- Reinforcing channels welded into larger models for additional strength to meet the highest expectations
- No open cavities to collect unwanted substances
- Cost-effective to meet your budget!

SPECIFICATIONS

- Dimensions: 33"L x 18"W x 34"H
- Construction: Stainless Steel, Solid Weld
- Casters: Ball Bearing
- Caster Diameter: 3"
- Warranty: Limited

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com