

# MCM 20" L x 4" W O.R. Table Width Extender | MCM 405



*Stock Photo - Options May Vary*

## **FACTORY NEW**

The MCM 405 O.R. Table Width Extender is designed to expand surgical table width by 4 inches, providing additional surface area and stability for wider surgical procedures. The unit measures 20 inches long by 4 inches wide and features integral Posi-Lock rail clamps and a sturdy aluminum side rail for secure attachment to existing table rails.

The extender locks securely to the operating table and allows quick setup without tools or complicated installation. Its intuitive design ensures consistent patient support during surgical procedures while maintaining workflow efficiency and safety standards required in the operating room.

Built for durability and convenience, the MCM 405 is ideal for facilities requiring expanded table capacity for specialized surgical cases. The pad (MCM 404) is sold separately.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)