

MCM Side Rail Extension | MCM 167

FACTORY NEW

The MCM-167 Side Rail Extension adds 6 inches of additional space to the end of operating room table rails, enabling easier placement of essential surgical accessories such as IV poles, anesthesia screens, and surgical lights.



Stock Photo - Options May Vary

Constructed from premium stainless steel, the extension provides lasting durability and corrosion resistance in demanding surgical environments. The solid construction ensures consistent performance throughout every procedure while maintaining easy cleaning and maintenance.

Installation is fast and tool-free, with the extension clamping directly to the end of the O.R. table rail without requiring permanent modifications. Compatible with U.S.-standard O.R. table rails, the MCM-167 works seamlessly with most surgical tables on the market.

FEATURES

- Adds 6 inches of additional rail space to the end of O.R. table
- Clamps securely to existing table rail with no tools required
- Constructed from high-quality stainless steel for durability
- Compatible with standard U.S. operating room table rails
- Supports various attachments, including IV poles, anesthesia screens, and surgical lights
- Enhances workflow by expanding accessory placement options
- Corrosion-resistant surface allows for easy cleaning and maintenance
- Cost-effective solution to increase table functionality without permanent modifications

SPECIFICATIONS

- Model: MCM-167
- Material: Stainless steel
- Extension Length: 6 inches added to existing rail
- Mounting Method: Clamps to end of OR table rail; no permanent modification required
- Compatibility: U.S.-standard OR table rails
- Origin: Made in USA

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)