

# MCM 25" Patient Transfer Roller Board | MCM 123



Stock Photo - Options May Vary

## FACTORY NEW

The MCM 123 Patient Transfer Roller Board is a 25-inch transfer device designed to assist healthcare staff with lateral patient transfers between surfaces. It reduces caregiver injury risk by eliminating the need for physical lifting during transfers between hospital beds, gurneys, operating room tables, X-ray tables, exam tables, and cast tables.

The board features smooth-operating rollers that enable effortless patient movement across surfaces. Its design prioritizes both functionality and ease of use in clinical and surgical environments. The device includes a cover and smooth-clean surface suitable for healthcare settings.

## FEATURES

- Designed to assist staff members with patient transfer
- Eliminates physical lifting of patients
- Allows effortless patient transfer from any two combinations of the following: hospital bed, gurney, operating room table, X-ray table, exam table, and cast table
- Smooth operating rollers make the transfer of patients easy
- Cover included.

## SPECIFICATIONS

- Length: 25"
- Width: 15"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)