


# MCM Patient Slider Board - 33"L | MCM 121



Stock Photo - Options May Vary

## FACTORY NEW

The MCM 121 patient slider board is manufactured from high-density polyethylene (HDPE) and measures 3/16" thick, offering an ideal balance of strength and flexibility. The board features a non-anti-static surface suitable for various clinical environments and includes a generous handle that provides a secure, ergonomic grip to simplify patient transfers.

Caregivers can easily transfer patients between a bed, gurney, stretcher, or OR table with this reliable board, reducing strain and improving safety. The hand holes double as openings for convenient storage and hanging. Designed for daily hospital use, the MCM 121 ensures smooth, efficient, and safe patient transfers.

## SPECIFICATIONS

- Length: 33"
- Width: 22"
- Height: 3/16"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)