

MCM Bariatric Patient Slider Board – 72”L x 32”W | MCM 122-W



Stock Photo - Options May Vary

FACTORY NEW

The MCM 122-W Bariatric Patient Slider Board reduces transfer-related risks while promoting safety for both staff and patients. Designed for smooth, stable movement, this board makes lateral transfers from stretcher to bed or operating table faster and safer.

Available in three sizes—half-length, full-length, and bariatric full-length—this board meets the needs of various care environments. Whether moving standard or bariatric patients, there is a size that fits each application.

Crafted from 3/16" thick HDPE, the board provides a durable, smooth-sliding surface that simplifies patient handling. Integrated hand holes make carrying and positioning easy, and allow staff to hang the board for storage, keeping it accessible and out of the way.

SPECIFICATIONS

- Length: 72"
- Width: 32"
- Height: Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com