


MCM Patient Slider Board - 72"L | MCM 122



Stock Photo - Options May Vary

FACTORY NEW

The MCM 122 is a premium-quality transfer board manufactured from high-density polyethylene (HDPE) measuring 3/16" thick. It provides a balance of strength and flexibility suitable for a wide range of clinical environments. The board features a generously sized handle that ensures a secure grip for efficient and safe patient transfers.

The MCM 122 slider board is ideal for transferring patients from bed to gurney, stretcher, or operating room table with minimal strain and maximum efficiency. Its smooth surface and robust construction help reduce the risk of patient injury or caregiver fatigue.

The board's hand holes serve a dual purpose, functioning as openings for convenient storage and hanging. The smooth-sliding HDPE surface facilitates easier transfers while maintaining durability for repeated clinical use.

SPECIFICATIONS

- Length: 72"
- Width: 22"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)