

MCM 30" Patient Transfer Roller Board | MCM 125



Stock Photo - Options May Vary

FACTORY NEW

The MCM 125 is a 30-inch patient transfer roller board designed to facilitate lateral patient transfers between surfaces such as beds, stretchers, and tables. It reduces physical strain on healthcare staff by eliminating the need to manually lift patients, improving safety for both caregivers and patients.

Equipped with smooth-operating rollers, the board allows effortless movement and can be used for transfers between any two combinations of hospital beds, gurneys, operating room tables, X-ray tables, exam tables, and cast tables. The board includes a protective cover and cleans easily for convenient maintenance.

Ideal for hospitals, clinics, and surgical centers, the MCM 125 supports efficient patient handling while prioritizing hygiene and safety in healthcare environments.

FEATURES

- Designed to assist staff members with patient transfer
- Eliminates physical lifting of patients
- Allows effortless patient transfer from any two combinations of the following: hospital bed, gurney, operating room table, X-ray table, exam table, and cast table
- Smooth operating rollers make the transfer of patients easy
- Cover included.

SPECIFICATIONS

- Length: 30"
- Width: 15"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)