

MCM 67" Patient Transfer Roller Board | MCM 130



Stock Photo - Options May Vary

FACTORY NEW

The MCM 130 Patient Transfer Roller Board measures 67 inches and is specifically designed to assist healthcare staff during lateral patient transfers. It eliminates the need for physical lifting, helping reduce the risk of caregiver injury while improving the overall safety of the transfer process. This roller board is ideal for use in hospitals, surgical centers, and long-term care facilities.

Smooth-operating rollers allow for effortless movement of patients between beds, stretchers, or tables, improving both comfort and efficiency. The board accommodates transfers between any two combinations of hospital beds, gurneys, operating room tables, X-ray tables, exam tables, and cast tables.

The MCM 130 includes a protective cover and features a smooth, easy-to-clean surface, ensuring hygiene standards are maintained in healthcare environments.

FEATURES

- Patient Transfer Roller Board: Full Body Length
- Designed to assist Staff Members with Patient Transfers
- Eliminates Physical Lifting of Patients
- Smooth Operating Rollers Makes the Transfer of Patients Easy
- Cover Included

SPECIFICATIONS

- Length: 67"
- Width: 15"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)