

MCM 30 x 16 Instrument Table with Shelf | MCM 502



Stock Photo - Options May Vary

FACTORY NEW

The MCM 502 Instrument Table is constructed from heavy-duty 16-gauge #304 stainless steel with solid welded construction throughout. The table features a 1.5" revealed edge, square tube legs with professionally welded and polished corners, and 3" ball-bearing casters for mobility.

Designed for durability and ease of cleaning, the MCM 502 includes no open cavities that could collect unwanted substances, making it ideal for sterile medical environments. The solid welded construction and square tube legs provide superior strength and longevity compared to alternative designs.

The table measures 30"L x 16"W x 34"H and comes with a limited 2-year warranty, providing a reliable instrument table solution for surgical and medical settings.

FEATURES

- Heavy Duty 16 gauge #304 stainless steel with a solid welded construction
- The revealed edge is 1-1/2" to give you the strongest table known
- Square tube legs allow more welding surface which has proven to outlast any round tube leg
- Corners are professionally welded and polished smooth
- 3" ball bearing casters are standard on all tables
- No open cavities to collect unwanted substances
- Reinforcing channels welded into larger models for additional strength to meet the highest expectations
- Cost-effective to meet your budget!

SPECIFICATIONS

- Dimensions: 30"L x 16"W x 34"H
- Construction: Stainless Steel, Solid Weld
- Casters: Ball Bearing
- Caster Diameter: 3"
- Warranty: Limited

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com