

MCM 33 x 18 Instrument Table w/Shelf | MCM 503



Stock Photo - Options May Vary

FACTORY NEW

The MCM 503 is a heavy-duty stainless steel instrument table with shelf, constructed from 16 gauge #304 stainless steel with solid welded construction. The table features a 1-1/2" revealed edge, square tube legs with professionally welded and polished corners, and 3" ball bearing casters for mobility. The design eliminates open cavities that could collect unwanted substances.

The square tube leg design provides more welding surface than round tube alternatives, resulting in superior durability and longevity. Reinforcing channels are welded into the structure for additional strength. All corners are professionally finished to ensure a smooth, hygienic surface suitable for medical environments.

This instrument table is designed to meet the demanding requirements of healthcare facilities while offering cost-effective value. The unit comes with a limited 2-year warranty.

FEATURES

- Heavy Duty 16 gauge #304 stainless steel with a solid welded construction
- The revealed edge is 1 1/2" to give you the strongest table known
- Square tube legs allow more welding surface which has proven to outlast any round tube leg
- Corners are professionally welded and polished smooth
- 3" ball bearing casters are standard on all tables
- Reinforcing channels welded into larger models for additional strength to meet the highest expectations
- No open cavities to collect unwanted substances
- Cost-effective to meet your budget!

SPECIFICATIONS

- Dimensions: 33"L x 18" x 34"H
- Construction: Stainless Steel, Solid Weld
- Casters: Ball Bearing
- Caster Diameter: 3"
- Warranty: Limited

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com