

MCM 36 x 20 Instrument Table w/Shelf | MCM 504



Stock Photo - Options May Vary

FACTORY NEW

The MCM 504 is a heavy-duty stainless steel instrument table with shelf, constructed from 16-gauge #304 stainless steel with solid welded construction. The table features a 1-1/2" revealed edge, square tube legs with professionally welded and polished corners, and 3" ball bearing casters for mobility. The design eliminates open cavities that could collect unwanted substances, making it suitable for clinical environments.

The table's square tube legs provide greater welding surface area compared to round tube designs, contributing to enhanced durability and longevity. All corners are professionally welded and polished smooth for a seamless finish. The solid welded construction ensures structural integrity and stability during use.

This instrument table is designed to meet the demands of healthcare facilities requiring reliable, easy-to-clean equipment. The stainless steel construction facilitates sanitation and resists corrosion, while the included shelf provides additional storage capacity for instruments and supplies.

FEATURES

- Heavy Duty 16 gauge #304 stainless steel with a solid welded construction
- The revealed edge is 1-1/2" to give you the strongest table known
- Square tube legs allow more welding surface which has proven to outlast any round tube leg
- Corners are professionally welded and polished smooth
- 3" ball bearing casters are standard on all tables
- Reinforcing channels welded into larger models for additional strength to meet the highest expectations
- No open cavities to collect unwanted substances
- Cost-effective to meet your budget!

SPECIFICATIONS

- Dimensions: 36"L x 20"W x 34"H
- Construction: Stainless Steel, Solid Weld
- Casters: Ball Bearing
- Caster Diameter: 3"
- Warranty: Limited

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com