

MCM 36 x 20 Instrument Table w/H Brace | MCM 513



Stock Photo - Options May Vary

FACTORY NEW

The MCM 513 is a heavy-duty instrument table constructed from 16-gauge #304 stainless steel with solid welded construction. It features a 1½" revealed edge, square tube legs, and professionally welded and polished smooth corners for maximum durability and strength.

The table is designed with no open cavities to collect unwanted substances, making it ideal for sterile medical environments. Standard 3" ball bearing casters provide smooth mobility, while the solid welded construction ensures long-lasting performance in demanding healthcare settings.

This instrument table combines robust construction with a practical design that meets the needs of operating rooms and surgical facilities. It comes with a limited 2-year warranty.

FEATURES

- Heavy Duty 16-gauge #304 stainless steel with a solid welded construction
- The revealed edge is 1-1/2" to give you the strongest table known
- Square tube legs allow more welding surface which has proven to outlast any round tube leg
- Corners are professionally welded and polished smooth
- 3" ball bearing casters are standard on all tables
- Reinforcing channels welded into larger models for additional strength to meet the highest expectations
- No open cavities to collect unwanted substances
- Cost-effective to meet your budget!

SPECIFICATIONS

- Dimensions: 36"L x 20"W x 34"H
- Construction: Stainless Steel, Solid Weld
- Casters: Ball Bearing
- Caster Diameter: 3"
- Warranty: Limited

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com