

MCM 48 x 20 Instrument Table w/Shelf | MCM 506



Stock Photo - Options May Vary

FACTORY NEW

The MCM 506 is a heavy-duty stainless steel instrument table with shelf, constructed from 16-gauge #304 stainless steel with solid welded construction. It features a 1½" revealed edge, square tube legs, and smooth polished corners for durability and professional appearance.

The table includes 3" ball bearing casters for mobility and is designed without open cavities to prevent buildup of unwanted substances. Square tube legs provide more welding surface than round alternatives, ensuring long-term reliability and structural integrity.

The MCM 506 is backed by a limited 2-year warranty and offers a cost-effective solution for healthcare facilities requiring reliable instrument storage and work surfaces.

FEATURES

- Heavy Duty 16 gauge #304 stainless steel with a solid welded construction
- The revealed edge is 1-1/2" to give you the strongest table known
- Square tube legs allow more welding surface which has proven to outlast any round tube leg
- Corners are professionally welded and polished smooth
- 3" ball bearing casters are standard on all tables
- Reinforcing channels welded into larger models for additional strength to meet the highest expectations
- No open cavities to collect unwanted substances
- Cost-effective to meet your budget!

SPECIFICATIONS

- Dimensions: 48"l x 20"w x 34"h
- Construction: Stainless Steel, Solid Weld
- Casters: Ball - Bearing
- Caster Diameter: 3"
- Warranty: Limited

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com