

# MCM 48 x 24 Instrument Table w/Shelf | MCM 507



Stock Photo - Options May Vary

## FACTORY NEW

The MCM 507 is a 48 x 24 inch instrument table with lower shelf storage, designed for surgical suites, hospitals, and outpatient facilities. Built from 16-gauge #304 stainless steel with solid welded construction, it provides a durable, spacious work surface that supports daily medical procedures.

The table features square tube legs that provide maximum welding surface for superior strength and longevity. All corners are professionally welded and polished smooth to eliminate sharp edges and ensure easy cleaning. The design includes no open cavities that could collect unwanted substances, maintaining the highest standards of sterility.

Mobility is enhanced by standard 3-inch ball bearing casters on all tables. The solid construction and precision engineering deliver the reliability and strength required in demanding healthcare environments.

## FEATURES

- Heavy Duty 16 gauge #304 stainless steel with a solid welded construction
- The revealed edge is 1-1/2" to give you the strongest table known
- Square tube legs allow more welding surface which has proven to outlast any round tube leg
- Corners are professionally welded and polished smooth
- 3" ball bearing casters are standard on all tables
- No open cavities to collect unwanted substances

## SPECIFICATIONS

- Dimensions: 48" L x 24" W x 34" H
- Construction: Stainless Steel, Solid Weld
- Casters: Ball Bearing
- Caster Diameter: 3"
- Warranty: Limited 2 Year

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)