

# MCM 36 x 24 Instrument Table w/H Brace | MCM 514



Stock Photo - Options May Vary

## FACTORY NEW

The MCM 514 is a heavy-duty stainless steel instrument table designed for surgical and clinical environments. Constructed from 16-gauge #304 stainless steel with solid welded construction, it features square tube legs with professionally welded and polished corners for maximum durability and strength. The table includes 3-inch ball bearing casters for mobility and a 1 1/2-inch revealed edge for structural integrity.

This instrument table is built without open cavities, eliminating areas where unwanted substances could collect and simplifying maintenance in demanding medical settings. The solid welded construction and reinforced design ensure reliable performance and long-term clinical use. The table measures 36 inches long by 24 inches wide by 34 inches high, making it suitable for a variety of surgical and instrument preparation applications.

The MCM 514 combines strength, mobility, and ease of maintenance in a cost-effective design. Its durable construction supports consistent performance in high-demand healthcare environments and comes with a limited 2-year warranty.

## FEATURES

- Heavy Duty 16 gauge #304 stainless steel with a solid welded construction
- The revealed edge is 1 1/2" to give you the strongest table known
- Square tube legs allow more welding surface which has proven to outlast any round tube leg
- Corners are professionally welded and polished smooth
- 3" ball bearing casters are standard on all tables
- Reinforcing channels welded into larger models for additional strength to meet the highest expectations
- No open cavities to collect unwanted substances

## SPECIFICATIONS

- Dimensions: 36"L x 24"W x 34"H
- Construction: Stainless Steel, Solid Weld
- Casters: Ball Bearing
- Caster Diameter: 3"
- Warranty: Limited

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)