

MCM 36 x 24 Instrument Table w/Shelf | MCM 505



Stock Photo - Options May Vary

FACTORY NEW

The MCM 505 is a 36 x 24 instrument table with shelf, constructed from 16-gauge #304 stainless steel with solid welded construction. Designed for hospitals, surgical centers, and clinical environments, it combines durability and functionality for reliable daily medical tasks.

The table features square tube legs that provide increased welding surface for superior strength and longevity compared to round tube designs. All corners are professionally welded and polished smooth to eliminate sharp edges and potential contamination points.

Equipped with 3-inch ball bearing casters for mobility and a bottom shelf for storage, the MCM 505 offers a no-cavity design that prevents unwanted substance accumulation. The solid construction and reinforced design ensure dependable support and long-lasting performance in demanding medical environments.

FEATURES

- Heavy Duty 16 gauge #304 stainless steel with a solid welded construction
- The revealed edge is 1 1/2" to give you the strongest table known
- Square tube legs allow more welding surface which has proven to outlast any round tube leg
- Corners are professionally welded and polished smooth
- 3" ball bearing casters are standard on all tables
- Reinforcing channels welded into larger models for additional strength to meet the highest expectations
- No open cavities to collect unwanted substances
- Cost-effective to meet your budget!

SPECIFICATIONS

- Dimensions: 24 x 36 x 34
- Construction: Stainless Steel, Solid Weld
- Casters: Ball Bearing
- Caster Diameter: 3"
- Warranty: Limited

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com