

MCM 48 x 24 Instrument Table w/H Brace | MCM 516



Stock Photo - Options May Vary

FACTORY NEW

MCM 516 Instrument Table w/H Brace. Mid Central Medical manufactures high-quality medical equipment to meet the demands of America's health care industry. MCM carries a Wide Variety of Innovative products, Including Stainless Steel instruments/Back Tables with Shelves in Many Sizes.

Additional MCM516 Features: This Heavy Duty Stainless Steel Instrument Table with H Brace is made from 16-gauge, #304 Stainless Steel with a Solid Welded Construction. Also featured, Revealed Edge is 1-1/2". Square Tube Legs, Corners Professionally Welded, and Polished Smooth. 3" Ball Bearing Casters are Standard on All of MCM's Tables. There are no open Cavities to Collect Unwanted Substances.

Moreover, the MCM 516 includes a limited 2-Year Warranty.

FEATURES

- Heavy Duty 16-gauge, #304 stainless steel with a solid welded construction
- The revealed edge is 1-1/2" to give you the strongest table known
- Square tube legs allow more welding surface which has proven to outlast any round tube leg
- Corners are professionally welded and polished smooth
- 3" ball bearing casters are standard on all tables

SPECIFICATIONS

- Dimensions: 48"l x 24"w x 34"h
- Construction: Stainless Steel, Solid Weld
- Casters: Ball Bearing
- Caster Diameter: 3"
- Warranty: Limited

**Brochure and images are a general overview of the product. Actual product options may vary.*