

# MCM 60 x 24 Instrument Table w/Shelf | MCM 508



*Stock Photo - Options May Vary*

## FACTORY NEW

The MCM 508 is a stainless steel instrument table designed for medical facilities requiring durable, cleanable work surfaces. Constructed from 16-gauge #304 stainless steel with solid welded construction, the table measures 60 inches long by 24 inches wide by 34 inches high and includes a bottom shelf for additional storage.

The table features square tube legs that provide maximum welding surface and proven durability compared to round tube designs. All corners are professionally welded and polished smooth to eliminate sharp edges and potential contamination points. The design includes no open cavities that could collect unwanted substances, making it ideal for sterile environments.

Standard 3-inch ball bearing casters enable easy mobility while the reinforced construction ensures stability and longevity in demanding healthcare settings. The table is backed by a limited 2-year warranty.

## FEATURES

- Heavy Duty 16-gauge #304 stainless steel with a solid welded construction
- The revealed edge is 1-1/2" to give you the strongest table known
- Square tube legs allow more welding surface which has proven to outlast any round tube leg
- Corners are professionally welded and polished smooth
- 3" ball bearing casters are standard on all tables
- Reinforcing channels welded into larger models for additional strength to meet the highest expectations
- No open cavities to collect unwanted substances

## SPECIFICATIONS

- Dimensions: 60" L x 24" W x 34" H
- Construction: Stainless Steel, Solid Weld
- Casters: Ball Bearing
- Caster Diameter: 3"
- Warranty: Limited

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)