

MCM 60 x 24 Instrument Table w/H Brace | MCM 517



Stock Photo - Options May Vary

FACTORY NEW

MCM 517 Instrument Table w/H Brace. Mid Central Medical manufactures high-quality medical equipment to meet the demands of America's health care industry. MCM carries a Wide Variety of Innovative products, Including Stainless Steel instruments/Back Tables with Shelves in Many Sizes.

Additional MCM517 Features: This Heavy Duty Stainless Steel Instrument Table with H Brace is made from 16-gauge, #304 Stainless Steel with a Solid Welded Construction. Also featured, Revealed Edge is 1-1/2". Square Tube Legs, Corners Professionally Welded, and Polished Smooth. 3" Ball Bearing Casters are Standard on All of MCM's Tables. Furthermore, there are no open Cavities to Collect Unwanted Substances—limited 2-Year Warranty.

FEATURES

- Heavy Duty 16-gauge #304 stainless steel with a solid welded construction
- The revealed edge is 1-1/2" to give you the strongest table known
- Square tube legs allow more welding surface which has proven to outlast any round tube leg
- Corners are professionally welded and polished smooth
- 3" ball bearing casters are standard on all tables

SPECIFICATIONS

- Construction: Stainless Steel, Solid Weld
- Table Dimensions: 24"wx 60"l x 34"h
- Casters: Ball Bearing
- Caster Diameter: 3"
- Warranty: Limited

**Brochure and images are a general overview of the product. Actual product options may vary.*