

# MCM 72 x 24 Instrument Table w/H Brace | MCM 518



Stock Photo - Options May Vary

## FACTORY NEW

The MCM 518 is a heavy-duty stainless steel instrument table designed for medical facilities. Constructed from 16-gauge #304 stainless steel with solid welded construction, it features square tube legs with professionally welded and polished corners for maximum durability and strength. The table includes a 1-1/2" revealed edge and is equipped with 3" ball bearing casters for mobility.

The design eliminates open cavities that could collect unwanted substances, making it ideal for sterile environments. The solid welded construction and reinforced design ensure reliable performance in demanding healthcare settings.

## FEATURES

- Heavy Duty 16 gauge #304 stainless steel with a solid welded construction
- The revealed edge is 1-1/2" to give you the strongest table known
- Square tube legs allow more welding surface which has proven to outlast any round tube leg
- Corners are professionally welded and polished smooth
- No open cavities to collect unwanted substances

## SPECIFICATIONS

- Table Dimensions: 24"w x 72"l x 34"h
- Construction: Stainless Steel, Solid Weld
- Casters: Ball Bearing
- Caster Diameter: 3"
- Warranty: Limited

*\*Brochure and images are a general overview of the product. Actual product options may vary.*