

MCM 72 x 24 Instrument Table w/Shelf | MCM 509



Stock Photo - Options May Vary

FACTORY NEW

MCM 509 Instrument Table w/Shelf. Mid Central Medical, an American company, manufactures high-quality medical equipment to meet the demands of America's healthcare industry. MCM carries a wide variety of innovative products, including stainless steel instrument/back tables with shelves in many sizes, designed for performance and reliability.

Additional MCM 509 Features: This heavy-duty table uses 16-gauge, #304 stainless steel with solid welded construction. It features a 1-1/2" revealed edge, square tube legs, smooth, polished corners, and 3" ball bearing casters.

MCM designs all tables with no open cavities to collect unwanted substances. Furthermore, the MCM 509 includes a limited 2-year warranty and is ideal for clinical, surgical, and procedural environments.

FEATURES

- Heavy Duty 16 gauge #304 stainless steel with a solid welded construction
- The revealed edge is 1-1/2" to give you the strongest table known
- Square tube legs allow more welding surface which has proven to outlast any round tube leg
- Corners are professionally welded and polished smooth
- 3" ball bearing casters are standard on all tables

SPECIFICATIONS

- Instrument Table Dimensions: 24"w x 72"l x 34"h
- Construction: Stainless Steel, Solid Weld
- Casters: Ball Bearing
- Caster Diameter: 3"
- Warranty: Limited

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)