

MCM Foot Operated Mayo Stand | MCM 751



Stock Photo - Options May Vary

FACTORY NEW

MCM Foot Operated Mayo Stand: MCM 751. Mid Central Medical manufactures high-quality medical equipment to meet the demanding standards of America's healthcare industry. As a result, the MCM 751 offers a strong balance of durability, functionality, and safety.

Moreover, this makes it a great choice for surgical centers, hospitals, and outpatient facilities. This stand features solid brass inserts that allow smooth, effortless height adjustment.

A built-in "Safe Descent" system ensures controlled lowering of the tray, protecting both patients and medical staff. Furthermore, the removable stainless steel tray adds convenience and simplifies cleaning between procedures.

FEATURES

- Low Profile Base Design
- Removable Tray
- Foot-Pedal Operated Lock
- Adjustable Height
- Easily Slides Under Equipment

SPECIFICATIONS

- Removable Tray Dimensions: 16" x 21"
- Adjustable Height Range: 36" to 62"
- Constructions: All Welded #304 Stainless Steel
- 2" Double Ball Bearing Casters
- 2 Year Limited Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)