

MCM Add-A-Rail | MCM 324



Stock Photo - Options May Vary

FACTORY NEW

The MCM 324 Add-A-Rail is specifically designed to fit hourglass-shaped surgical tables and provides versatility and reliable performance in clinical environments. Constructed from solid stainless steel, it offers long-lasting strength and resistance to corrosion in busy surgical settings.

The MCM 324 mounts securely to blocks underneath the table using two countersunk screws, creating a stable fit while keeping the table surface clear and free from obstructions. This easy installation helps staff work more efficiently during procedures.

The open-end design allows accessories to slide on quickly from the end, and users can attach multiple accessories simultaneously to adjust for different procedures. This flexibility enables staff to work faster and more safely during surgery.

FEATURES

- For Hourglass-shaped tables only
- Solid Stainless Steel
- Mounts to blocks underneath the table
- Easy Mounting - Two Countersunk Screws
- Allows Accessories to be slipped on from the end
- Allows You to Adapt Multiple Accessories at one time

SPECIFICATIONS

- For Hour Glass Shaped Tables
- Dimensions: 24"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)