

MCM Hourglass Arm and Hand Surgery Table | MCM 320



Stock Photo - Options May Vary

FACTORY NEW

The MCM 320 Hourglass Arm/Hand Surgery Table is engineered for arm and hand surgical procedures with a distinctive hourglass design that expands to 24 inches at the ends, providing surgeons with improved ergonomic access and comfort during operations. The table features a single adjustable post leg with a 2-inch-thick pad and attaches securely to standard operating tables using exclusive Posi-Lock Clamps.

Constructed with a radiolucent surface, the MCM 320 delivers outstanding C-Arm imaging quality without compromising stability or positioning. The lightweight design combines strength with portability, while the hourglass shape accommodates a wide variety of arm and hand surgical procedures. A built-in handle enables effortless carrying and serves as a convenient storage bracket.

FEATURES

- Phenolic Hourglass Arm and Hand Surgery Table
- Single Post Leg
- Unique Posi-Lock Clamps
- Comes with 2" Thick Pad
- Hourglass Shape Expands the Table Width to 24" at the ends to Enhance the Surgeon's Comfort.
- Adjustable Leg
- A built-in handle allows effortless carrying
- Radiolucent offering exceptional C-arm imaging

SPECIFICATIONS

- Weight: 10 lbs.
- Operating Surface Dimension: 15" x 32"
- Table Width: Expands to 24" at the ends
- The leg adjustment has a range of 12"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)