

MCM Hourglass Arm and Hand Surgery Table | MCM 321



Stock Photo - Options May Vary

FACTORY NEW

The MCM 321 is a phenolic hourglass arm and hand surgery table designed for a wide range of surgical procedures. Its distinctive hourglass shape expands to 24 inches at the ends to enhance surgeon comfort and access. The table features a 2-inch thick pad and secure Posi-Lock Clamps for stable attachment to the operating table.

The MCM 321 combines lightweight construction with exceptional durability and stability. Its radiolucent design provides excellent C-arm imaging quality, making it ideal for procedures requiring fluoroscopic guidance. The table accommodates various surgical approaches and patient positioning needs.

Key design elements include a Tee Foot for positioning, an adjustable leg with 12 inches of range, and a built-in handle for easy transport and convenient storage. The operating surface measures 15 by 32 inches and weighs only 10 pounds, making it practical for busy surgical environments.

FEATURES

- Phenolic Arm/Hand Table
- Tee Foot
- Unique Posi-Lock Clamps
- Comes with 2" Thick Pad
- Operating surface Accommodate a Wide Variety of Procedures
- Adjustable Leg
- The built-in handle allows effortless carrying
- Radiolucent offering exceptional C-arm imaging

SPECIFICATIONS

- Weight: 10 lbs.
- Operating Surface Dimension: 15" x 32"
- Table Width: Expands to 24" at the ends
- The leg adjustment has a range of 12"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)