

MCM Hourglass Arm and Hand Surgery Table | MCM 326



Stock Photo - Options May Vary

FACTORY NEW

The MCM 326 is a carbon fiber arm and hand surgery table designed for secure positioning during surgical procedures. It features Posi-Lock clamps that attach firmly to the operating room table, providing unmatched stability. The extra-thick carbon fiber construction delivers superior strength while maintaining a lightweight design at just 12 pounds.

The operating surface measures 15" x 32" and accommodates a wide variety of procedures. The distinctive hourglass shape expands to 24" width at the ends, enhancing surgeon comfort and access. A 2" thick pad provides cushioning for patient positioning.

The table is fully radiolucent for clear imaging during procedures. An adjustable leg and built-in handle enable effortless carrying and convenient storage. Optional Add-A-Rail accessories can be attached to expand functionality for multiple surgical applications.

FEATURES

- Carbon Fiber Table
- Clamps positively to operating table with our unique Posi-lock clamps
- Made with extra-thick carbon fiber to provide superior strength when it is most essential
- Totally radiolucent for extremely clear images
- The operating surface is 15" x 32" to accommodate a wide variety of procedures
- Hourglass shape expands the table width to 24" at the ends to enhance the surgeon's comfort
- Includes 2" thick pad
- The built-in handle allows effortless carrying and also acts as a convenient storage bracket, providing a place to hang the unit from when in storage

SPECIFICATIONS

- Weight: 12 lbs.
- Operating surface: 15" x 32"
- Table Width: 24" at the ends
- Leg adjustment Range: 12"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)