


MCM Add-A-Rail | MCM 314



Stock Photo - Options May Vary

FACTORY NEW

The MCM 314 Add-A-Rail is a stainless steel accessory rail designed for rectangular-shaped surgical tables. It mounts securely to blocks underneath the table using two countersunk screws, allowing accessories to be slipped on from the end without additional tools or complex installation.

Constructed from solid stainless steel, the MCM 314 provides a durable, corrosion-resistant solution that resists wear and rust. Its minimalist design integrates seamlessly with existing table setups while maintaining a polished, modern appearance.

The rail enables simultaneous adaptation of multiple accessories, maximizing table functionality and flexibility. The straightforward mounting process and user-friendly design make it an efficient solution for optimizing surgical table capabilities.

FEATURES

- For Rectangular-shaped table only
- Solid Stainless Steel
- Mounts to blocks underneath the table
- Easy Mounting - Two Countersunk Screws
- Allows Accessories to be slipped on from the end

SPECIFICATIONS

- For Rectangle Shaped Tables
- Dimensions: 18"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)