

MCM Arm and Hand Surgery Table | MCM 310



Stock Photo - Options May Vary

FACTORY NEW

The MCM 310 is an arm and hand surgery table designed for surgical precision and patient safety. It features a single post leg with a comfortable 2" thick pad for patient support, adjustable leg support for versatility, and reinforced aluminum bars along the outer edges for added strength and stability during procedures.

The table secures firmly to the operating table using integrated Posi-Lock Clamps, ensuring minimal movement during use. A built-in handle allows effortless carrying and serves as a convenient storage bracket when not in use.

The MCM 310 is radiolucent, allowing excellent C-Arm imaging without interference. Built for strength, comfort, and precision, it accommodates a wide range of arm and hand procedures in any surgical setting.

FEATURES

- Single Post Leg
- Comes complete with 2" Pad
- Adjustable Leg
- Integrated Unique Posi-Lock Clamps
- Reinforced Outside Edges provide added Strength
- Operating Surface accommodates a Wide Variety of Procedures
- Built-in Handle allows Effortless Carrying
- Radiolucent, offering Excellent C-Arm Imaging

SPECIFICATIONS

- Operating Surface: 18" x 32"
- Leg Adjustment Range: 12"
- Aluminum Bar Reinforcement: 3/4" X 3/4" (at the very outside of table)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)