

MCM Rectangle Arm / Hand Surgery Table | MCM 311



Stock Photo - Options May Vary

FACTORY NEW

The MCM 311 is a rectangular arm and hand surgery table designed for strength, precision, and flexibility across a variety of surgical procedures. It features a sturdy Tee Foot base and includes a 2-inch-thick comfort pad for enhanced patient support. The generous operating surface accommodates hand surgeries and minor arm operations with stable, ergonomic positioning.

The table secures firmly to the operating table using MCM's exclusive Posi-Lock Clamps, which prevent movement during procedures and ensure maximum safety. An adjustable leg provides customization to meet specific surgical requirements. Aluminum bars reinforce the outer edges for added durability and strength during heavy-duty clinical use.

The MCM 311 is radiolucent and fully compatible with C-arm imaging for clear, unobstructed views. A built-in handle allows effortless carrying and serves as a convenient storage bracket, combining portability with practical functionality.

FEATURES

- Phenolic Arm/Hand Table
- Tee Foot
- Adjustable Leg
- Clamps Positively to Operating Table
- Comes complete with 2" Thick Pad
- Integrated Unique Posi-Lock Clamps
- Reinforced Outside Edges provide added Strength
- Radiolucent, offering Excellent C-Arm Imaging

SPECIFICATIONS

- Operating Surface: 18" x 32"
- Leg Adjustment Range: 12"
- Aluminum Bar Reinforcement: 3/4" x 3/4" (at the very outside edges)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)