

# MCM Adjustable Anesthesia Screen | MCM 101



*Stock Photo - Options May Vary*

## **FACTORY NEW**

The MCM 101 Adjustable Anesthesia Screen is constructed from high-quality stainless steel for durability and easy maintenance in operating room environments. It features a 5/8" diameter mounting post that securely fits into a standard Clark socket (sold separately as MCM 160), enabling quick setup and stable positioning on most surgical tables.

The internal clamp adjusts easily from 16" to 26" in height, allowing users to set the screen at the appropriate level for each procedure. This adjustable range improves patient privacy while maintaining visibility for the anesthesia provider during surgical operations.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)