

MCM Uni-Pole IV Pole | MCM 216



Stock Photo - Options May Vary

FACTORY NEW

The MCM 216 Uni-Pole IV Pole is a universal clamping IV pole designed to secure to operating table rails, table legs, bed rails, and similar surfaces. It clamps to any round or square tube from 7/8" to 1-1/4" in diameter, making it adaptable to most healthcare environments.

Constructed of 304 stainless steel and aluminum with composite handles, the pole resists rust and corrosion even in harsh operating room conditions and is easy to clean. The main clamp secures the pole to furniture, while a second handle enables 360-degree rotation for flexible positioning.

The pole's height adjusts from 18" to 36" via a third handle. The sturdy 2-hook top features ram horn-style hooks that securely hold IV bags without risk of slipping. All components are integrated into one unit, eliminating the need for separate clamps or attachments.

FEATURES

- Stainless Steel construction
- Unique clamp mounts to most surfaces up to 1 1/4" diameter
- Adjusts from 18" to 36"
- Extremely universal
- No need for additional clamps

SPECIFICATIONS

- 2 prong hooks
- Hooks: Ram Horn Style
- Clamp to any round or square tube from 7/8" to 1-1/4"
- Height Adjustments: 18" to 36" in height

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)