

# MCM Free Access Arm and Hand Surgery Table | MCM 301



Stock Photo - Options May Vary

## FACTORY NEW

The MCM 301 is a free access arm and hand surgery table designed to provide surgeons with optimal visibility and easy access to the surgical field. Its unique shape eliminates leg obstruction, allowing unrestricted C-arm and X-ray procedures. The radiolucent design supports state-of-the-art imaging technologies.

The operating surface measures 12" x 30" and is suitable for fistula, angiogram, and angioplasty procedures. Unique Posi-lock clamps securely attach to the operating room table, providing a positive and stable hold during procedures.

The table ships complete with a pad and clamps, offering a ready-to-use solution for surgical teams. Its durable construction and economical design make it a practical addition to any operating room.

## FEATURES

- Allows easy access
- No leg to restrict C-arm or X-ray procedures
- Radiolucent - allowing state-of-the-art C-arm access
- Clamps positively to operating table with unique Posi-lock clamps
- The operating surface is 12" x 30"
- Great for many procedures, including Fistula, Angiogram, and Angioplasty
- Economical to fit your budget
- Extremely durable

## SPECIFICATIONS

- Weight: 10 lbs
- Width: 12"
- Length: 30"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)