

MCM Carpal Arm and Hand Surgery Table | MCM 330



Stock Photo - Options May Vary

FACTORY NEW

The MCM 330 is a carpal arm and hand surgery table featuring a legless design that provides excellent rigidity while allowing easy access for C-arm imaging. The 18" x 32" operating surface is one of the largest available, yet the table weighs only 11 pounds, offering a spacious work area without compromising portability or stability.

The MCM 330 clamps securely to the OR rail using Posi-Lock Clamps, providing extra stability during delicate procedures. It is compatible with any American-made operating table and ships complete with a 2" pad and clamps for immediate surgical use.

Ideal for carpal tunnel surgeries, the MCM 330 is versatile enough for fistula, angiogram, and angioplasty procedures. Its combination of size, weight, and secure positioning makes it a reliable choice for hand and arm surgical applications.

FEATURES

- Legless design
- Optimum rigidity
- Easy access for C-arm procedures
- One of the largest Operating Surfaces available
- Clamps positively to operating table rail with unique Posi-lock clamps for added stability
- The ultimate choice for carpal tunnel surgery, among other procedures
- Compatible with any American made operating table
- Table Ships Complete with Pad

SPECIFICATIONS

- Weight: 11 lbs
- Width: 18"
- Length: 32"
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)