

MCM Carry Cart | MCM 145



Stock Photo - Options May Vary

FACTORY NEW

MCM Carry Cart: MCM 145. Mid Central Medical, a trusted American manufacturer, produces high-quality medical equipment to support healthcare professionals across the country. Known for durability and innovation, MCM designs each product to meet real-world clinical demands.

The MCM 145 is a practical tool that supports daily efficiency and cleanliness in medical settings. The MCM 145 is specifically designed to work with the MCM 140 Interlocking Step Stool, but also accommodates similar equipment.

Additionally, it allows staff to easily lift, transport, and store step stools, reducing strain and keeping floors clear for safer, easier cleaning. Its compact frame helps save space while keeping equipment organized and accessible.

FEATURES

- Stainless Steel
- Fully Welded Construction
- Provides Storage/Transportation for Stacking Stools
- 3" Ball Bearing Neoprene Casters
- Designed to Carry 4 - 6 Stools

SPECIFICATIONS

- Construction: 16 Gauge Stainless Steel
- Casters: 3" Scuff-Resistant Casters
- Dimensions: 14"D x 21"W x 39"H
- Limited 2 Year Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)