

# MCM Central Sterile Instrument Wrapping Station | MCM 596



Stock Photo - Options May Vary

## FACTORY NEW

MCM Central Sterile Instrument Wrapping Station: MCM 596. Mid Central Medical delivers dependable solutions for sterile processing with the MCM 596. Built with a stainless steel top and upper rack, it provides a strong, stable surface for wrapping instruments in clean, controlled environments.

The ergonomic 35" work height reduces strain on staff, making long tasks more comfortable and efficient in high-demand healthcare settings. This wrapping station features heavy-duty 4" locking casters that allow smooth movement and secure positioning.

Workers can easily relocate the unit and lock it firmly in place when needed. Its adjustable size fits a variety of workspace layouts, offering the flexibility essential for central sterile departments and other sterile processing areas.

## FEATURES

- Sturdy work style table design.
- Removable bins for easy cleaning and reloading instruments.
- Easy to move but can still lock in place with locking casters.

## SPECIFICATIONS

- Construction: Stainless Steel
- Table Dimensions: 32"D x 60"L x 34"H
- Work Height: 35"
- 4" Heavy-Duty, Locking Caster

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)