

MCM Chrome 5-Leg 2 Hook IV Pole | MCM 201



Stock Photo - Options May Vary

FACTORY NEW

The MCM 201 is a chrome IV pole with a sturdy 5-leg base designed for stable fluid management in hospitals, clinics, and surgical centers. The 2-hook top securely holds IV bags while the 24-inch diameter base with rubber bumpers provides enhanced stability during patient transport or stationary use. The pole features a 1-inch diameter chrome construction with height adjustment from 52 to 93 inches.

This IV pole is compatible with a range of attachments for customization to specific medical procedures and patient care requirements. The 2-inch high scuff-resistant casters enable smooth mobility in inpatient units, outpatient settings, and surgical suites. Accessories are sold separately, allowing users to adapt the pole to fit changing clinical needs.

The MCM 201 delivers reliable performance and dependable operation backed by a limited 2-year warranty. For facilities requiring additional hook capacity, a 4-hook version (MCM 202) is also available.

FEATURES

- 2 Hook Top
- 24" Diameter 5-Leg Epoxy Coated Base
- 2" High Scuff-Resistant Casters
- 1" Diameter Chrome Pole with Friction-Lock Knob
- Height Adjustable from 52"-93"
- Rubber Bumpers on End of Base Legs
- Very Durable and Dependable
- Limited 2 Year Warranty

SPECIFICATIONS

- 2 Hook Top
- Base Diameter: 24"
- Caster Size: 2" High
- Pole Diameter: 1"
- Height: 52"-93"
- Limited 2 Year Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*