

MCM Chrome 5-Leg 4 Hook IV Pole | MCM 202



Stock Photo - Options May Vary

FACTORY NEW

The MCM 202 Chrome IV Pole features a 5-leg base and 4-hook top designed for stable, reliable fluid management in hospitals and clinics. Its durable chrome finish and wide base provide maximum stability even when fully loaded with multiple IV bags.

The pole is height adjustable from 52 to 93 inches and includes a 1-inch diameter chrome pole with friction-lock knob for secure positioning. The 24-inch diameter base is epoxy coated with rubber bumpers on the end of each leg to protect flooring.

The MCM 202 supports multiple IV bags simultaneously, making it ideal for patients requiring complex or multiple fluid treatments. A 2-hook version is also available as the MCM 201, and IV pole attachments are sold separately for customization.

FEATURES

- 24" Diameter 5-Leg Epoxy Coated Base
- 2" High Scuff-Resistant Casters
- 1" Diameter Chrome Pole with Friction-Lock Knob
- Height Adjustable from 52"-93"
- Rubber Bumpers on End of Base Legs
- Limited 2 Year Warranty
- Also Available with 2-Hook Top
- Attachments Sold Separately

SPECIFICATIONS

- 4 Hook Top
- Base Diameter: 24"
- Caster Size: 2" High
- Pole Diameter: 1"
- Height: 52"-93"
- Limited 2 Year Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*