

# MCM Chrome 5-Leg/4 Hook I.V. Pole | MCM 273



Stock Photo - Options May Vary

## FACTORY NEW

The MCM 273 is a chrome IV pole featuring a 5-leg spider base with a 4-hook top designed for secure fluid administration in clinical settings. The adjustable pole extends from 50 to 92 inches in height and is controlled by a friction knob for easy positioning. Scuff-resistant casters provide smooth mobility while the chrome construction offers durability and easy cleaning.

The pole comes standard with a charcoal base and is available in five base colors upon request: yellow, green, blue, charcoal, and red. A 2-hook top version is also available to accommodate different departmental needs. The unit includes a 2-year limited warranty.

IV pole attachments are sold separately to customize the setup for specific clinical applications. The MCM 273 is ideal for hospitals, clinics, and other healthcare facilities requiring reliable fluid administration support.

## FEATURES

- 5-Leg Spider Leg 3"
- 4 Hook Top
- Scuff Resistant Casters for Ease of Mobility
- Chrome Pole
- Adjustable Pole with Friction Knob
- Standard Base Color is Charcoal
- Base Available in 5 Colors Upon Request
- Limited 2 Year Warranty

## SPECIFICATIONS

- Base Weight: 22 Pounds
- Caster Size: 3"
- Chrome Pole Diameter: 1"
- Height: 50"-92"
- Warranty: Limited 2 Year

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)